










 <p>WATER CHAMBER 1 every 6 months</p>	 <p>FULL FACE MASK (Covers nose & mouth) Mask w/ headgear: 1 every 6 months Mask w/o headgear: 1 every 3 months</p>	 <p>NASAL MASK (Covers Nose) Mask w/ headgear: 1 every 6 months Mask w/o headgear: 1 every 3 months</p>	 <p>NASAL PILLOWS Mask w/ headgear: 1 every 6 months Mask w/o headgear: 1 every 3 months</p>
 <p>DISPOSABLE FILTERS 2 per month</p>	 <p>FULL FACE MASK SEAL/CUSHION 1 per month</p>	 <p>NASAL MASK SEAL/CUSHION 2 per month</p>	 <p>NASAL PILLOWS 2 per month</p>
 <p>REUSABLE FILTER 1 every 6 months</p>	 <p>HEADGEAR 1 every 6 months</p>	 <p>HEATED TUBING 1 every 3 months</p>	 <p>STANDARD TUBING 1 every 3 months</p>