



SAFETY PRECAUTIONS FOR OXYGEN CONCENTRATOR USE

Please note:

Oxygen is not combustible or flammable on its own. Oxygen "supports combustion", which means that a fire will ignite more readily and burn more rapidly in the presence of oxygen.

Safety Guidelines:

1. Do not use or store oxygen by an open flame.
2. Do not under ANY circumstances smoke while wearing your oxygen.
3. Do not cook over a gas stove with your oxygen on.
4. Do not lay concentrator flat when transporting, always keep upright.
5. Do not use extension cords on concentrator.
6. If there are pets in the home, check tubing regularly to make sure there are no punctures or bite marks.
7. If there are area rugs on slippery floors, please remove.
8. Keep equipment in a well-ventilated area at least six (6) inches from wall.
9. Always post "NO SMOKING" sign in area where oxygen is used.
10. We strongly recommend the use of smoke detectors.
11. Use the oxygen only at prescribed liter flow and hours.
12. Keep oxygen equipment on flat surface.